

# Chapter 10

## Simple Enzyme Experiments

*Sharon B. Miller*

Science Department  
Tri-County Technical College  
Pendleton, South Carolina 29670

Sharon Miller received her B.A. (1967) from Winthrop College in biology, M.S. (1969) from Clemson University in plant pathology, and Ph.D. (1974) from the University of Georgia in botany/mycology. She has taught biology, botany, and/or mycology for 15 years at either the high school or college level, and served 2 years as laboratory coordinator for the multi-section introductory biology course for non-majors at Clemson University. At present, she is an adjunct professor at Clemson University and a part-time biology instructor at Tri-County Technical College. Her special interests include laboratory instruction and working with public school teachers. She was one of five Project WILD facilitators in South Carolina to receive a meritorious service award in 1986.

**Reprinted from:** Miller, S. B. 1992. Simple enzyme experiments. Pages 153-161, *in* Tested studies for laboratory teaching, Volume 6 (C.A. Goldman, S.E. Andrews, P.L. Hauta, and R. Ketchum, Editors). Proceedings of the 6<sup>th</sup> Workshop/Conference of the Association for Biology Laboratory Education (ABLE), 161 pages.

- Copyright policy: <http://www.zoo.utoronto.ca/able/volumes/copyright.htm>

Although the laboratory exercises in ABLE proceedings volumes have been tested and due consideration has been given to safety, individuals performing these exercises must assume all responsibility for risk. The Association for Biology Laboratory Education (ABLE) disclaims any liability with regards to safety in connection with the use of the exercises in its proceedings volumes.

## Contents

Introduction.....	154
Experiments	
Amylase .....	154
Catalase.....	155
Catecholase .....	156
Invertase.....	157
Papain.....	158
Pectinase .....	158
Pepsin .....	159
Rennin.....	160
Literature Cited.....	161

## Introduction

Thousands of enzymes are found in living cells where they act as catalysts for the thousands of chemical reactions which occur. In addition to making life possible, many enzymes have numerous applications that affect our daily lives in other ways such as food processing, clinical diagnoses, sewage treatment, and the textile industry. Enzyme experiments are ideal for “hands on” opportunities and since several factors affect the rate at which enzymatic reactions proceed, an enzyme experiment presents many opportunities in the biology laboratory. The experiments presented in this chapter were selected on the basis that they: (1) are simple to prepare, (2) use materials which are familiar to the student, and (3) can be used as a base from which to construct additional experiments. Instructions are given for a basic reaction which can then be modified to investigate various aspects of enzyme activity such as the effect of temperature, enzyme concentration and substrate concentration, depending on the level at which they are used. The following enzymes are included: amylase, catalase, catecholase, invertase, papain, pectinase, pepsin, and rennin.

Except for the pepsin experiment, all experiments can be completed during a 2- to 3-hour laboratory period. The time for each individual experiment varies from “instant” results with catalase to an incubation period of 30 minutes or more with amylase. Except for the enzyme solutions, materials listed with each experiment are based on what one student will need if working alone. The enzyme solution is usually enough for a class. Most quantities are expressed in more than one way so that the experiments can be done with a minimum of equipment. For further reading see Allan et al. (1983), Mader (1982), Morholt et al. (1966), Novo Laboratories (1975), and Porter et al. (1973).

## Amylase

### *Introduction*

Amylase is an enzyme that catalyses the hydrolysis of the polysaccharide starch (amylose) to the disaccharide maltose. It is readily abundant in saliva, but somewhat unpleasant to obtain in large quantities. It is widely distributed in plant tissues, but is most abundant in seeds, where it apparently functions in initiating the breakdown of stored starch to glucose which is needed in large amounts during germination.

*Materials* starch-agar plates (0.2% soluble starch, 2% agar)

Wax pencil  
Distilled water (in wash bottle)

### *Procedure*

1. Prepare starch-agar plates (do not have to be sterile if used within a day or two). Allow to solidify and cool.
2. Use a wax pencil to label the bottom of the plate: “soaked seeds”, “boiled seeds”, “dry seeds”, etc. (You might want to include a few drops of saliva from your mouth for comparison.)
3. Use a sharp razor blade to cut the corn grain longitudinally and place, cut surface down, onto the agar surface. (You may wish to dissect out the embryo.) Be sure to space corn grains at least 2 cm from each other.
4. Incubate for 30 minutes.
5. Remove corn and rinse plate gently with distilled water.
6. Flood plate with iodine solution, swish around as color develops, rinse with distilled water, record results. (Any clear areas of agar can be removed and tested for sugars.)

### *Results*

After flooding the plates with iodine solution, the agar will stain a deep purple in all areas where starch remains. Areas of agar where dead seeds were placed will be purple, likewise for dry seeds (unless the incubation period is much longer) since dormant seeds produce very little amylase. Areas of the agar covered by saliva, or by a living embryo, will appear clear since the starch has been broken down.

## **Catalase**

### *Introduction*

Hydrogen peroxide ( $\text{H}_2\text{O}_2$ ) is naturally formed in living organisms, however it is very harmful and is broken down immediately by several enzymes including catalase. This enzyme catalyses the breakdown of hydrogen peroxide to water and oxygen. Persons with acatalasemia (a hereditary condition) have extremely low catalase activity and, although present worldwide, it is more commonly found in Koreans.

### *Materials*

Test tubes (one for each material to be tested plus extra for control)  
Hydrogen peroxide ( $\text{H}_2\text{O}_2$ ) (3% solution)  
Assorted living tissue: sliced raw potato, ground meat, liver, yeast cells, ground young leaves  
Assorted non-living material: piece of baked potato or cooked liver, etc. (Use caution with rocks or sand, some will “bubble.”)

### *Procedure*

1. Fill each labelled test tube approximately 1/3 full with fresh hydrogen peroxide.
2. Add a small amount of material to be tested.
3. Note whether or not bubbles are produced.

### *Results*

Any fresh, living material will normally have enough catalase present to produce bubbles of gas (oxygen) upon exposure to the hydrogen peroxide. Dry yeast and liver are most impressive when used in this experiment.

## **Catecholase**

### *Introduction*

The brown color, which usually develops when a potato or apple is cut or bruised, is the result of a chemical catalyzed by several enzymes. One of the most important of these enzymes is catecholase (catechol oxidase or polyphenol oxidase). This enzyme is brought into contact with its substrate, catechol, when cells are ruptured and exposure to air (oxygen) resulting in the oxidation of catechol to the brownish colored benzoquinone. This so called “wound reaction” is apparently protective since the quinones are toxic to microorganisms.

### *Materials*

0.1% Catechol solution (in dropper bottle with dark glass)  
Distilled water (150 ml plus full dropper bottle)  
Electric blender  
Test tubes (2)  
Test tube rack  
Dropper bottle for potato extract  
Beaker (250 ml)  
Cheese cloth  
Stirring rods (2)  
Balance

### *Procedure*

1. Prepare catecholase by thoroughly blending potatoes (15 g/100 ml of water) in distilled water. Filter through cheese cloth. Put filtrate in dropper bottle.
2. Fill each of two test tubes 1/4 full (3 ml) of distilled water.
3. Add 10 drops of catechol to each tube.
4. Add 10 drops of potato extract to one tube, label C.
5. Add 10 drops of distilled water to the unlabelled test tube. Mix.
6. Note and record color of each tube at 1-minute intervals for 5–10 minutes.

### *Results*

The contents of tube C (catechol plus enzyme) will darken with time as the reaction proceeds while the control tube remains clear.

## **Invertase**

### *Introduction*

Our most common food sugar—the disaccharide, sucrose—is formed in all green plants. The metabolism of sucrose in the animal body begins with the action of invertase (sucrase) which hydrolyzes the disaccharide to two monosaccharides, fructose and glucose. This same enzyme is also produced by plants and fungi.

### *Materials*

Sucrose (0.25 M solution)  
Glucose (0.25 M solution)  
Benedict's solution (in dropper bottle)  
Distilled water (100 ml)  
Hot plate with water bath  
Stirring rods (2)  
Beaker (50 ml)  
Test tubes (5)  
Test tube rack  
Pipet or syringe to measure 3 ml  
Tes-Tape (tape from drug store used to test for glucose in urine samples)  
Balance or teaspoon

### *Procedure*

1. Prepare yeast by mixing 1 teaspoon (3 g) dry yeast with 20 ml of distilled water. Let stand for 20 minutes.
2. Fill each of two test tubes 1/3 full with sucrose solution.
3. Add 3 ml yeast suspension to one tube (label 1). Mix.
4. Add 3 ml distilled water to the other tube. Mix.
5. After 10 minutes test each test tube plus the yeast suspension with a strip of Tes-Tape. (Benedict's test for reducing sugars can also be used here since sucrose will give a negative Benedict's test and glucose/fructose will give a positive test (yellow-orange-red) depending on the amount of reducing sugar present. Remove about 2 ml of the solution, place in another test tube, add 10 drops of Benedict's solution and place in a boiling water bath for about 3 minutes. Prepare a glucose solution for comparison.)

### *Results*

For best results, do not go over the recommended incubation period. Both the Tes-Tape and Benedict's give positive tests for samples from tube 1 while the tube without yeast gives negative results.

## Papain

### *Introduction*

Papain, from the latex of the papaya plant, is one of a family of plant enzymes that includes bromelin (from pineapple) and ficin (from fig), all of which break down proteins. This is why the directions on a box of Jello remind you never to use fresh or frozen pineapple in your gelatin, since gelatin is the protein responsible for the “gel.” A convenient source of papain is fresh pineapple juice or meat tenderizer.

### *Materials*

Gelatin (Knox, Jello)  
Beaker (150 ml)  
Balance or teaspoon  
Stirring rods (3)  
Test tubes (2)  
Test tube rack  
Beaker of ice water  
Hot plate  
Distilled water (100 ml)

### *Procedure*

1. Prepare a gelatin solution by heating 1 teaspoon (3.0 g) of gelatin in 100 ml distilled water until dissolved. (Gently mix, do not boil.) Cool to room temperature.
2. Pour meat tenderizer into one of the two test tubes until it fills approximately 0.5 cm of the tube. Label this tube as P. Do not put meat tenderizer in the other tube.
3. Fill each test tube 1/3 full (5 ml) with the gelatin solution. Mix gently.
4. Place tubes in ice water for 10 minutes.
5. Remove from ice bath and note the degree of gelatinization.

### *Results*

The tube without meat tenderizer (papain) will contain firm gelatin. Tube P which contains papain will be almost liquid.

## Pectinase

### *Introduction*

Pectinases are enzymes that breakdown the polysaccharide pectin which is located primarily in the middle lamella of cell walls. Pectin compounds are widely distributed in plant tissues, especially in fruits. They are large, colloidal molecules which are responsible for holding dispersed particles in suspension in fruit juices. In addition to influencing the amount of suspended particles in a juice, pectins also increase the viscosity of the juice. The presence of suspended particles in tomato and orange juice is not undesirable to most people, however, “clear” juices such as apple and grape are more often preferred. Commercially prepared pectinase can be added to prepared fruits in order to hasten the release of juice and aid in the “settling out” of suspended particles in fruit juice.

*Materials*

Funnels (2)  
 Graduated cylinders, 100-ml (2)  
 Beakers, 50-ml (2)  
 Pipets or syringe, to measure 0.5 ml liquid (2)  
 Spatulas or spoons (2)  
 Apple sauce  
 Distilled water (10 ml)  
 Pectinase (available from several biological supply companies or businesses providing cider and/or wine making supplies)  
 Wax pencil

*Procedure*

Cheesecloth (two squares to fit funnel)

1. Place approximately 25 ml of apple sauce into each of the two beakers labelled “no enzyme” and “pectinase”.
2. Add 0.5 ml distilled water to the “no enzyme” beaker and 0.5 ml pectinase to the “pectinase” beaker.
3. Stir the contents of each beaker thoroughly (using separate spatulas).
4. Let stand 10 minutes or longer.
5. Place cheesecloth in a funnel and place the funnel into the graduated cylinder.
6. With the aid of the spatula, pour contents of each beaker into a separate funnel and collect filtrate.
7. Record the amount of juice collected in each cylinder.

*Results*

The volume of filtrate from the apple sauce plus pectinase is usually at least double that of the filtrate without the enzyme.

**Pepsin***Introduction*

Pepsin is an enzyme which aids in the breakdown of proteins and is secreted in the stomach of most animals. It functions at a very low pH, which makes it ideal to use for studying the effect of pH on enzyme activity. Commercial or natural albumin can be used as the protein source; however hard-boiled egg white is more impressive to demonstrate the action of pepsin.

*Materials*

Egg white (hard-boiled)  
 Test tubes (6)  
 Test tube rack

Solutions:

- A: 1.0% Pepsin
- B: Pepsin (1.0%) in 0.4% hydrochloric acid
- C: 0.4% Hydrochloric acid
- D: Pepsin (1.0%) in 0.5% sodium bicarbonate
- E: 0.5% Sodium bicarbonate
- F: Distilled water

Metric ruler

Knife

*Procedure*

1. Label test tubes A to F.
2. Fill each tube 1/3 full (5 ml) with the corresponding solution.
3. Drop a small (2 mm) cube of egg white into each tube.
4. Incubate at room temperature for approximately 12 hours. (The speed of this reaction can be increased by using very thin strips of egg white and/or incubating at 30°C. Students might want to compare the action of papain or bromelin under similar conditions.)
5. Examine tubes for the presence of the egg white.

*Results*

After 12 hours, the egg white will be “digested” in tube B; little, if any, of the cube will remain, while there appears to be little change in the other tubes. If allowed to incubate longer, the cube in tube D and possibly C will decrease in size.

## Rennin

*Introduction*

Rennin is an enzyme obtained from the fourth stomach of a calf or other young bovines. It causes rapid clotting of milk by causing certain bonds to break in the soluble casein molecule (milk protein) converting it to the insoluble casein, thus producing “curdled” milk. Most other proteolytic enzymes will curdle milk, but they also cause continued breakdown of the casein, unlike rennin. Other forms of rennin exist as rennet (an extract of the enzymes from the fourth stomach, i.e., not pure rennin) and rennilase or hannilase which are trade names for milk-clotting microbial enzymes. This enzyme is available at some health food stores and is present in Junket which is found with the home-made ice cream supplies at many grocery stores.

*Materials*

Junket

Milk (do not use UHT long shel-life type)

Distilled Water (50ml)

Beaker (100 ml)

Stirring rods (3)

Test tubes (2)

Test tube rack

3. To one of the tubes add an equal volume of the Junket solution (label R). Mix.
4. Add an equal volume of distilled water to the other tube. Mix.
5. If the milk is cold, let tubes sit for 30 minutes or more before recording results. If the milk has been brought to room temperature before steps 3 and 4, then examine tubes after 10 minutes.

### Result

The tube (R) containing the enzyme will start to clot while the other tube will remain liquid.

### Literature Cited

- Allan, M. A., S. Aneja, F. R. Armstrong, R. A. Garcia, D. R. Helms, R. N. Holmes, S. B. Miller, A. D. Smith, D. J. Stroup, W. M. Surver, and C. K. Wagner. 1983. *General biology: A laboratory manual for Biology 105*. Contemporary Publishing Co., Raleigh, North Carolina, 265 pages. (See pages 95–104 for a detailed study of catecholase.)
- Mader, S. S. 1982. *Laboratory manual for Inquiry Into Life*. Third edition. W. C. Brown Co., Dubuque, Iowa, 335 pages. (Pages 38–41 provide experiments with catalase, urease, and rennin.)
- Morholt, E., P. F. Brandwein, and A. Joseph. 1966. *A sourcebook for the biological sciences*. Second edition. Harcourt, Brace, and World Inc., New York, 795 pages. (Pages 206–217 include experiments using amylase, invertase, papain, pepsin, rennin, and others.)
- Novo Laboratories. 1975. *Enzymes: Nature's catalyst. Teacher's Manual*. Novo Laboratories, Wilton, Connecticut, 40 pages. (Available from Carolina Biological Supply Co., 2700 York Rd., Burlington, NC 27215. An excellent, although brief, source of basic information and numerous simple experiments.)
- Porter, D., W. M. Darley, E. L. Dunn, A. J. Hicks, and H. W. Rines. 1973. *Botany laboratory manual*. Paje Publishing, Atlanta, Georgia. (Pages 17–26 include experiments using amylase and phosphorylase from plants.)